

Recent Items Navigator Plugin

USER MANUAL

User manual for Recent Items Navigator 1.3

Windows, Mac OS X and Linux

February 15, 2017

This software is for research purposes only.

QIAGEN Aarhus
Silkeborgvej 2
Prismet
DK-8000 Aarhus C
Denmark



Contents

1	Introduction to the Recent Items Navigator	4
2	Install	5
3	Uninstall	7

Chapter 1

Introduction to the Recent Items Navigator

Recent Items Navigator can be used to see the elements that you have recently opened or used.

To see the recent items, click the tab to the left of the **Navigation Area** (🌐) (see figure 1.1). Double-clicking the item will open it in a view.




Figure 1.1: Recent items.

When you restart, the list of recent items will be reset and built from scratch.

There can be max 50 items in the list - when this limit has been reached, the oldest item will be removed when a new is added.

Chapter 2

Install

Recent Items Navigator is installed as a plugin. Plugins are installed using the plugin manager. In order to install plugins on Windows, the Workbench must be run in administrator mode: Right-click the program shortcut and choose "Run as Administrator". Then follow the procedure described below.

Help in the Menu Bar | Plugins... ()

or **Plugins () in the Toolbar**

The plugin manager has two tabs at the top:

- **Manage Plugins.** This is an overview of plugins that are installed.
- **Download Plugins.** This is an overview of available plugins on CLC bio's server.

To install a plugin, click the **Download Plugins** tab. This will display an overview of the plugins that are available for download and installation (see figure 2.1).

Clicking a plugin will display additional information at the right side of the dialog. This will also display a button: **Download and Install.**

Click the Recent Items Navigator and press **Download and Install.** A dialog displaying progress is now shown, and the plugin is downloaded and installed.

If the Recent Items Navigator is not shown on the server, and you have it on your computer (for example if you have downloaded it from our website), you can install it by clicking the **Install from File** button at the bottom of the dialog. This will open a dialog where you can browse for the plugin. The plugin file should be a file of the type ".cpa".

When you close the dialog, you will be asked whether you wish to restart the CLC Workbench. The plugin will not be ready for use until you have restarted.




Figure 2.1: The plugins that are available for download.

Chapter 3

Uninstall

Plugins are uninstalled using the plugin manager:

Help in the Menu Bar | Plugins... ()

or **Plugins () in the Toolbar**

This will open the dialog shown in figure 3.1.




Figure 3.1: *The plugin manager with plugins installed.*

The installed plugins are shown in this dialog. To uninstall:

Click the Recent Items Navigator | Uninstall

If you do not wish to completely uninstall the plugin but you don't want it to be used next time you start the Workbench, click the **Disable** button.

When you close the dialog, you will be asked whether you wish to restart the workbench. The plugin will not be uninstalled until the workbench is restarted.