

# **Recent Items Navigator Plugin**

USER MANUAL

# User manual for Recent Items Navigator 1.4

Windows, Mac OS X and Linux

November 7, 2017

**This software is for research purposes only.**

QIAGEN Aarhus  
Silkeborgvej 2  
Prismet  
DK-8000 Aarhus C  
Denmark



# Contents

|          |   |          |
|----------|---|----------|
| <b>1</b> | <b>Introduction to Recent Items Navigator</b> | <b>4</b> |
| <b>2</b> | <b>Install and uninstall plugins</b>          | <b>5</b> |
| 2.1      | Install . . . . .                             | 5        |
| 2.2      | Uninstall . . . . .                           | 6        |

# Chapter 1

## Introduction to Recent Items Navigator

Recent Items Navigator can be used to see the elements that you have recently opened or used.

To see the recent items, click the tab to the left of the **Navigation Area** (🌐) (see figure 1.1). Double-clicking the item will open it in a view.




Figure 1.1: Recent items.

When you restart, the list of recent items will be reset and built from scratch.

There can be max 50 items in the list - when this limit has been reached, the oldest item will be removed when a new is added.



## Chapter 2

# Install and uninstall plugins

Recent Items Navigator is installed as a plugin.

**Note:** In order to install plugins and modules, the Workbench must be run in administrator mode. On Linux and Mac, it means you must be logged in as an administrator. On Windows, you can do this by right-clicking the program shortcut and choosing "Run as Administrator".

Plugins are installed and uninstalled using the plugin manager.

**Help in the Menu Bar | Plugins... (  ) or Plugins (  ) in the Toolbar**

The plugin manager has two tabs at the top:

- **Manage Plugins.** This is an overview of plugins that are installed.
- **Download Plugins.** This is an overview of available plugins on QIAGEN Aarhus server.

### 2.1 Install

To install a plugin, click the **Download Plugins** tab. This will display an overview of the plugins that are available for download and installation (see figure 2.1).




Figure 2.1: The plugins that are available for download.

Select Recent Items Navigator to display additional information about the plugin on the right side

of the dialog. Click **Download and Install** to add the plugin functionalities to your workbench.

### Accepting the license agreement

Part of the installation involves checking and accepting the end user license agreement (EULA) as seen in figure 2.2.




Figure 2.2: Read the license agreement carefully.



Please read the EULA text carefully before clicking in the box next to the text **I accept these terms** to accept. If requested, fill in your personal information before clicking **Finish**.

If Recent Items Navigator is not shown on the server but you have the installer file on your computer (for example if you have downloaded it from our website), you can install the plugin by clicking the **Install from File** button at the bottom of the dialog and specifying the plugin \*.cpa file saved on your computer.

When you close the dialog, you will be asked whether you wish to restart the workbench. The plugin will not be ready for use until you have restarted.

## 2.2 Uninstall

Plugins are uninstalled using the plugin manager:

**Help in the Menu Bar** | **Plugins...** (  ) or **Plugins** (  ) in the **Toolbar**

This will open the dialog shown in figure 2.3.

The installed plugins are shown in the **Manage plugins** tab of the plugin manager. To uninstall, select Recent Items Navigator and click **Uninstall**.

If you do not wish to completely uninstall the plugin, but you do not want it to be used next time you start the Workbench, click the **Disable** button.




Figure 2.3: *The plugin manager with plugins installed.*

When you close the dialog, you will be asked whether you wish to restart the workbench. The plugin will not be uninstalled until the workbench is restarted.